


Product Name: Recombinant Human ZCCHC7 Protein
Catalog #: PEH10329



Summary

Name	ZCHC7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N3Z6
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NSESSSKPN QKKLIVLSDS EVIQLSDGSE VITLSEDSI YRCKGKNVRV QAQENAHGLS SSLQSNELVD KKCKSDIEKP KSEERSGVIR EVMIIIEVSSS EEEEESTISEG DNVESWMLLG CEVDDKDDDI LLNLVGCENS

image



Background

Alternative Names	Zinc finger CCHC domain-containing protein 7 (TRAMP-like complex RNA-binding factor ZCCHC7)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use