

Product Name: Recombinant Human ARID3B Protein
Catalog #: PEH10330



Summary

Name	ARI3B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IVW6
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QQQQKQPHL APLQMDAREK QGQQMREAQF LYAQLVTQP TLLSATAGRP SGSTPLGPLA RVPPTAAVAQ VFERGNMNSE PEEEDGGLED EDGDDEVAEV AEKETQAASK YFHVQKVARQ

image



Background

Alternative Names AT-rich interactive domain-containing protein 3B (ARID domain-containing protein 3B) (Bright and dead ringer protein) (Bright-like protein)

Note

For Research Use Only , Not for Diagnostic Use