


Product Name: Recombinant Human FBH1 Protein
Catalog #: PEH10331



Summary

Name	FBH1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NFZ0
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RWTNRDPNHG LYPKPRTKRG SRGQGSQRCI PEFFLAGKQP CTNDMAKSNS VGQDSCQDSE GDMIFPAESS CALPQEGSAG PGSPGSAPPS RKRSWSSEEE SNQATGTSRW DGVSKKAPRH HLSVPCTRPR EARQEAEDST SRLSAESGET

image



Background

Alternative Names F-box DNA helicase 1 (hFBH1) (EC 3.6.4.12) (F-box only protein 18)

Note

For Research Use Only , Not for Diagnostic Use