


Product Name: Recombinant Human ZHX1 Protein
Catalog #: PEH10334



Summary

Name	ZHX1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UKY1
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TRAESSIDE EVHESVSDSN QQNKKVEGGY ECKYCTFQTP DLNMFTFHVD SEHPNVVLNS SYVCVECNFL TKRYDALSEH NLKYHPGEEN FKLTMVKRNN QTIFEQTIND LTFDGSFVKE ENAEQAESTE VSSSGISISK TPIMKMMKKN VENKRIAVHH NSVEDVPEEK ENEIKPDREE IVENPSSAS

image



Background

Alternative Names Zinc fingers and homeoboxes protein 1

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZHX1 Protein
Catalog #: PEH10334

