


Product Name: Recombinant Human WDR12 Protein
Catalog #: PEH10338



Summary

Name	WDR12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9GZL7
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DNKKYAVDDV PFSIPAASEI ADLSNIINKL LKDKNEFHKK VEFDFLIKQQ FLRMPLDKHM EMENISSEEV VEIEYVEKYT APQPEQCMFH DDWISSIKGA EEWILTGSYD KTSRIWSLEG KSIMTIVGHT DVVKDVAWVK KDSLSCLLS ASMDQTILLW EWNVERNKVK ALHCCRGHAG SVDSIAVDGS GTKFCSGSWD KMLKIWSTVP TDEEDEMEES TNRPKKQKT

image



Background

Alternative Names Ribosome biogenesis protein WDR12 (WD repeat-containing protein 12)

Note

Product Name: Recombinant Human WDR12 Protein
Catalog #: PEH10338



For Research Use Only , Not for Diagnostic Use