


Product Name: Recombinant Human RAPGEF5 Protein
Catalog #: PEH10339



Summary

Name	RPGF5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92565
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PHLERKDSAA ALSDRELPLP TFDVPYFKYI DEEDEDDEWS SRSQSSTEDD SVDSLSDRY VVSGTPEKI LEHLLNDLHL EEVQDKETET LLDDFLTYT VFMTTDDLQ ALLRHYSACK YQGKEENSDV PRRKRKVLHL VSQWIALYKD WLPEDEHSM FLKTIYRNVL DDVYEYPILE KELKEFQKIL

image



Background

Alternative Names	Rap guanine nucleotide exchange factor 5 (Guanine nucleotide exchange factor for Rap1) (M-Ras-regulated Rap GEF) (MR-GEF) (Related to Epac) (Repac)
--------------------------	---

Note

Product Name: Recombinant Human RAPGEF5 Protein
Catalog #: PEH10339



For Research Use Only , Not for Diagnostic Use