


Product Name: Recombinant Human SYCP2 Protein
Catalog #: PEH10340



Summary

Name	SYCP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BX26
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DNIDKHIKTA KCVENTENKN VEFNQNFSE LQDVIPDSQA AEKRDHTILP GVLDNICGNK IHSKWACWTP VTNIELCNNQ RASTSSGDTL NQDIVINKKL TKQKSSSSIS DHNSEGTGKV KYKKEQTDHI KIDKAEVEVC

image



Background

Alternative Names	Synaptonemal complex protein 2 (SCP-2) (Synaptonemal complex lateral element protein) (hsSCP2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use