


Product Name: Recombinant Human ZIC5 Protein
Catalog #: PEH10343



Summary

Name	ZIC5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96T25
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HKRTHTGEKP FKCEFDGCDR KFASSDRKK HSHVHTSDKP YYCKIRGCDK SYTHPSSLRK HMKIHCKSPP PSPGPLGYSS VGTPVGAPLS PVLDPARSHS STLSPQVTNL NEWYVCQASG APSHLHTPSS NGTTSETEDE EIYGNPEVVR

image



Background

Alternative Names Zinc finger protein ZIC 5 (Zinc finger protein of the cerebellum 5)

Note

For Research Use Only , Not for Diagnostic Use