


Product Name: Recombinant Human IP6K2 Protein
Catalog #: PEH10344



Summary

Name	IP6K2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UHH9
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PEVVLDSDAE DLEDLSEESA DESAGAYAYK PIGASSVDVR MIDFAHTTCR LYGEDTVVHE GQDAGYIFGL QSLIDIVTEI

image



Background

Alternative Names Inositol hexakisphosphate kinase 2 (InsP6 kinase 2) (InsP6K2) (EC 2.7.4.-)
(P(i)-uptake stimulator) (PiUS)

Note

For Research Use Only , Not for Diagnostic Use