


Product Name: Recombinant Human BRF1 Protein
Catalog #: PEH10347



Summary

Name	TF3B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92994
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LGISDSIREC ISSQSSDPKD ASGDGELDLS GIDDLEIDRY ILNESEARVK AELWMRENAE YLREQREKEA RIAKEKELGI YKEHKPKKSC KRREPIQAST AREAIEKMLE QKKISSKINY

image



Background

Alternative Names	Transcription factor IIIB 90 kDa subunit (TFIIIB90) (hTFIIIB90) (B-related factor 1) (BRF-1) (hBRF) (TAF3B2) (TATA box-binding protein-associated factor, RNA polymerase III, subunit 2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human BRF1 Protein
Catalog #: PEH10347

