


Product Name: Recombinant Human CCNB1IP1 Protein
Catalog #: PEH10355



Summary

Name	CIP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NPC3
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GYAWVTACSH IFCDQHSGE FSRSPAICPA CNSTLSGKLD IVRTELPSE EYKAMVLAGL RPEIVLDISS RALAFWYQV HQERLYQEYN FSKAEGHLKQ MEKIYTQQIQ SKDVELTSMK GEVTSMKKVL EEYKKKFSDI SEKLMERNRQ YQKLQGLYDS LRLRNITIAN HEGTLEPSMI AQSGVLGFPL GNNSKFPLDN TPVRNRGDGD GDFQFRPFFA GSPTAPEPSN SFFSFVSPSR

image



Background

Alternative Names	E3 ubiquitin-protein ligase CCNB1IP1 (EC 2.3.2.27) (Cyclin-B1-interacting protein 1) (Human enhancer of invasion 10) (RING-type E3 ubiquitin transferase CCNB1IP1)
--------------------------	--

Product Name: Recombinant Human CCNB1IP1 Protein
Catalog #: PEH10355



Note

For Research Use Only , Not for Diagnostic Use