


Product Name: Recombinant Human ZFP91 Protein
Catalog #: PEH10362



Summary

Name	ZFP91_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96JP5
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RHRDTENTRS SRSKTGSLQL ICKSEPNTDQ LDYDVGEEHQ SPGGISSEEE EEEEEEMLIS EEEIPFKDDP RDETYKPHLE RETPKPRRKS GKVKEEKEKK EIKVEVEVEV

image



SDS



WB

Background

Alternative Names	E3 ubiquitin-protein ligase ZFP91 (EC 2.3.2.27) (RING-type E3 ubiquitin transferase ZFP91) (Zinc finger protein 757) (Zinc finger protein 91 homolog) (Zfp-91)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZFP91 Protein
Catalog #: PEH10362

