


Product Name: Recombinant Human DDX23 Protein
Catalog #: PEH10363



Summary

Name	DDX23_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BUQ8
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EHGDKKPKAQ PLSLEELLAK KKAEEEEAEAK PKFLSKAERE AEALKRRQQE VEERQRMLEE ERKKRKQFQD LGRKMLEDPO ERERRERRER MERETNGNED EEGRQKIREE KDKSKELHAI KERYLGGIKK RRRTRHLNDR KVFVEWDASE DTSIDYNPLY KERHQVQLLG

image



Background

Alternative Names	Probable ATP-dependent RNA helicase DDX23 (EC 3.6.4.13) (100 kDa U5 snRNP-specific protein) (DEAD box protein 23) (PRP28 homolog) (U5-100kD)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human DDX23 Protein
Catalog #: PEH10363

