


Product Name: Recombinant Human ZFH4 Protein
Catalog #: PEH10366



Summary

Name	ZFH4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86UP3
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AVAEDDEKDT SERDNSEGKN SNKDSVSVAG GTQPLLLAKE EDVATKRSKP TEDNKFCEHQ FYQCPYCNYN SRDQSRIQMH VLSQHSVQPV ICCPLCQDVL SNKMHLQLHL THLHSVSPDC VEKLLMTVPV PDVMMPNSML LPAAASEKSE RDTAAAVTAE GSGKYSGESP MDDKSMAGLE DSKANVEVKV EEQKPTKEPL EVSEWNKNSS KDVKIPDTLQ

image



Background

Alternative Names Zinc finger homeobox protein 4 (Zinc finger homeodomain protein 4) (ZFH-4)

Note

Product Name: Recombinant Human ZFHX4 Protein
Catalog #: PEH10366



For Research Use Only , Not for Diagnostic Use