



Product Name: Recombinant Human PHF21A Protein
Catalog #: PEH10367



Summary

Name	PF21A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96BD5
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PGTRLAGPQT VQLSKPSLEK QTVKSHTETD EKQTESRTIT PPAAPKPKRE ENPQKLAFMV SLGLVTHDHL EEIQSKRQER KRRTTANPVY SGAVFEPERK KSAVTYLNST MHPGTRKRGR PPKYNAVLGF GALTPTSPQS SHPDSPENEK TETTFTFPAP VQPVSLPSPT STDGDIHEDF

image



Background

Alternative Names PHD finger protein 21A (BHC80a) (BRAFF35-HDAC complex protein BHC80)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PHF21A Protein
Catalog #: PEH10367

