


Product Name: Recombinant Human KCTD13 Protein
Catalog #: PEH10369



Summary

Name	BACD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WZ19
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AAAPSLEAPK PSGLEPGPAA YGLKPLTPNS KYVKLNVGGS LHYTTLRTLT GQDTMLKAMF SGRVEVLTDG GGWVLIDRSG RHFGTILNYL RDGSVPLPES TRELGELLGE ARYYLVQGLI EDCQLALQKQ RETLSPLCLI PMVTSPREEQ QLLASTSKPV VKLLHNRSNN KYSYTSTSDG NLLKNIELFD KLALRFHGRL LFLKDVLGDE ICCWSFYGQG RKIAEVCCTS IVYATEKKQT KVEFPEARIF EETLNILIYE TPRGPDPALE EATGGAAGAG GAGRGEDEEN

image



Background

Alternative Names	BTB/POZ domain-containing adapter for CUL3-mediated RhoA degradation protein 1 (hBACURD1) (BTB/POZ domain-containing protein KCTD13) (Polymerase delta-interacting protein 1) (TNFAIP1-like protein)
--------------------------	--

Product Name: Recombinant Human KCTD13 Protein
Catalog #: PEH10369



Note

For Research Use Only , Not for Diagnostic Use