

Product Name: Recombinant Human PRDM13 Protein
Catalog #: PEH10375



Summary

Name	PRD13_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H4Q3
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AAAALPGARY AQLPPAPGLP LERCALPLD PGGLKAYPGG ECSHLPVAMP AFTVYNGELL YGSPATTAYY PLKLHFGGLL KYPESISYFS GPAAAALSPA ELGSLASIDR EIAMHNQQLS EMAAGKGRGR LDSGTLPPAV AAAGGTGGGG SGGSGAGKPK

image



Background

Alternative Names PR domain zinc finger protein 13 (EC 2.1.1.-) (PR domain-containing protein 13)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PRDM13 Protein
Catalog #: PEH10375

