


Product Name: Recombinant Human TAF5 Protein
Catalog #: PEH10378



Summary

Name	TAF5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15542
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LQEKQNNQIW NIVQEHLIYD IFDGMPSKQ QIDAMVGSLA GEAKREANKS KVFFGLLKEP EIEVPLDDED EGENEEGKP KKKKPKKDSI GSKSKKQDPN APPQNRIPLP ELKDSDKLDK IMNMKETTKR VRLGPDCLPS ICFYTFLNAY QGLTAVDVTD DSSLIAGGFA DSTVRVWSVT PKKLRSVKQA

image



Background

Alternative Names	Transcription initiation factor TFIID subunit 5 (Transcription initiation factor TFIID 100 kDa subunit) (TAF(II)100) (TAFII-100) (TAFII100)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TAF5 Protein
Catalog #: PEH10378

