


Product Name: Recombinant Human CHD6 Protein
Catalog #: PEH10383



Summary

Name	CHD6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TD26
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RKSETTVESL ELDQGLTNPS LRSPEESTES TDSQKRRSGR QVKRRKYNED LDFKVVDDDG ETIAVLGAGR TSALSASTLA WQAEPPEDD ANIEKILAS KTVQEVHPGE PPFDELEFYV KYRNFSYLHC KWATMEELEK DPRIAQKIKR FRNKQAQMKH IFTEPDEDLF NPDYVEVDRI LEVAHTKDAE TGEEVTHYLV KWCSLPYEEES TWELEEDVDP AKVKEFESLQ VLPEIKHVER PASDSWQKLE KSREYKNSNQ

image



Background

Alternative Names	Chromodomain-helicase-DNA-binding protein 6 (CHD-6) (EC 3.6.4.12) (ATP-dependent helicase CHD6) (Radiation-induced gene B protein)
--------------------------	--

Product Name: Recombinant Human CHD6 Protein
Catalog #: PEH10383



Note

For Research Use Only , Not for Diagnostic Use