


Product Name: Recombinant Human SNRK Protein
Catalog #: PEH10384



Summary

Name	SNRK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NRH2
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PARAADSVLN GHRSKGLCDS AKKDDLPELA GPALSTVPPA SLKPTASGRK CLFRVEEDEE EDEEDKKPMS LSTQVVLRRK PSVTNRLTSR KSAPVLNQIF EEGESDDEFD MDENLPPKLS RLKMNIASPG TVHKRYHRRK

image



Background

Alternative Names SNF-related serine/threonine-protein kinase (EC 2.7.11.1) (SNF1-related kinase)

Note

For Research Use Only , Not for Diagnostic Use