


Product Name: Recombinant Human LDB2 Protein
Catalog #: PEH10385



Summary

Name	LDB2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43679
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EPMQELMSRH KTYNLSPDC LKTCLEQKWQ RMVAPPAEPT RQPTTKRRKR KNSTSSTSNS SAGNNANSTG SKKKTAAANL SLSSQVPDVM VVGEPTLMGG EFGDEDERLI TRLENTQYDA ANGMDDEEDF

image



Background

Alternative Names	LIM domain-binding protein 2 (LDB-2) (Carboxyl-terminal LIM domain-binding protein 1) (CLIM-1) (LIM domain-binding factor CLIM1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use