


Product Name: Recombinant Human FAAP24 Protein
Catalog #: PEH10386



Summary

Name	FAP24_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BTP7
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PVHVPLGHIV ANEKWRGSQL AQEMQGKIKL IFEDGLTPDF YLSNRCCILY VTEADLVAGN GYRKLVRVR NSNNLKGIVV VEKTRMSEQY FPALQKFTVL DLGMVLLPVA SQMEASCLVI QLVQEQTKEP

image



Background

Alternative Names	Fanconi anemia core complex-associated protein 24 (Fanconi anemia-associated protein of 24 kDa)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use