


**Product Name: Recombinant Human RNASEH2C Protein**  
**Catalog #: PEH10388**



## Summary

<b>Name</b>	RNH2C_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q8TDP1
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	19 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	MESGDEAAIE RHRVHLRSAT LRDAVPATLH LLPCEVAVDG PAPVGRFFTP AIRQGPEGLE VSRGRCLRG EEVAVPPGLV GYVMVTEEEK VSMGKPDPLR DSGTDDQEEE PLERDFDRFI GATANFSRFT LWGLETIPGP DAKVRGALTW PSLAAAIHAQ VPED

## image



## Background

<b>Alternative Names</b>	Ribonuclease H2 subunit C (RNase H2 subunit C) (Aicardi-Goutieres syndrome 3 protein) (AGS3) (RNase H1 small subunit) (Ribonuclease H1 subunit C)
--------------------------	---

## Note

**Product Name: Recombinant Human RNASEH2C Protein**  
**Catalog #: PEH10388**



---

For Research Use Only , Not for Diagnostic Use