



Product Name: Recombinant Human DPF3 Protein
Catalog #: PEH10389



Summary

Name	DPF3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92784
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KKRRLHPPED PKLRLLEIKP EVELPLKKG FTSESTTLEA LLRGEGVEKK VDAREEESIQ EIQRVLENDE NVEEGNEED LEEDIPKRKN RTRGRARGSA GRRRRHDAAS QEDHDKPYVC DICGKRYKNR PGLSYHYAHT HLAEEGDEA QDQETRSPPN

image



Background

Alternative Names	Zinc finger protein DPF3 (BRG1-associated factor 45C) (BAF45C) (Zinc finger protein cer-d4)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human DPF3 Protein
Catalog #: PEH10389

