

Product Name: Recombinant Human MIS18BP1 Protein
Catalog #: PEH10391



Summary

| | |
|------------------------------|---|
| Name | M18BP_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q6P0N0 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 30 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | SSLKLNHKK NQFLKMTTFN NKNIFQSTML TEATTSNSSL DISAIKPNKD GLKNKANYES PGKIFLRMKE KVL RDKQEQP SRNSSLLEPQ KSGNNETFTP NRVEKKKLQH TYLCEEKENN KSFQDDSSL RASVQGVPLE SSNNDIFLPV KQKIQCQQEK KAPLHNLTYE LPTLNQEQEN FLAVEARNKT LTRAQLAKQI FHSKESIVAT TSKKDTFVL ESVDSADEQF QNTNAETLST NCIPKNGSL |

image



Background

Alternative Names Mis18-binding protein 1 (Kinetochores-associated protein KNL-2 homolog) (HsKNL-2) (P243)

Note

Product Name: Recombinant Human MIS18BP1 Protein
Catalog #: PEH10391



For Research Use Only , Not for Diagnostic Use