



Product Name: Recombinant Human RNF113A Protein
Catalog #: PEH10392



Summary

Name	R113A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15541
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RDSGKQKAAY GDSSEEEEE NEPESLGVVY KSTRSAKPVG PEDMGATAVY ELDTEKERDA QAIFERSQKI QEELRGKEDD KIYRGINNYQ KYMKPKDTSM GNASSGMVRK GPIRAPEHLR ATVRWDYQPD ICKDYKETGF CGFGDSCKFL HDRSDYKHGW QIERELDEGR YGVYEDENYE VGSDDEEIPF

image



Background

Alternative Names	E3 ubiquitin-protein ligase RNF113A (EC 2.3.2.27) (Cwc24 homolog) (RING finger protein 113A) (Zinc finger protein 183)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RNF113A Protein
Catalog #: PEH10392

