


Product Name: Recombinant Human HELZ2 Protein
Catalog #: PEH10393



Summary

Name	HELZ2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BYK8
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GDEALSPASR DITATTAQTE AAAAPAGDAV KEDVVPGACA AGAAAAAGVE STEAEADAEAD FWPWDGELNA DDAILRELLD ESQKVMVTVG EDGLLDTVAR PESLQQARLY ENLPPAALRK LLHAEPERYR HCSFVPETFE RASAIPLDDA SSGPIQVRGR LDCGMAFAGD EVLVQLLSGD KAPEGRLRGR

image



Background

Alternative Names	Helicase with zinc finger domain 2 (ATP-dependent helicase PRIC285) (Helicase with zinc finger 2, transcriptional coactivator) (PPAR-alpha-interacting complex protein 285) (PPAR-gamma DNA-binding domain-interacting protein 1) (PDIP1) (PPAR-gamma DBD-interacting protein 1) (Peroxisomal proliferator-activated receptor A-interacting complex 285 kDa)
--------------------------	--

Product Name: Recombinant Human HELZ2 Protein
Catalog #: PEH10393



protein) (EC 3.6.4.-)

Note

For Research Use Only , Not for Diagnostic Use