


Product Name: Recombinant Human DIP2A Protein
Catalog #: PEH10395



Summary

Name	DIP2A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14689
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	APLPAEVRES LAELELESE GDITQKGYEK KRAKLLARYI PLIQGIDPSL QAENRIPGPS QTAAAPKQQ KSRPTASRDE RFRSDVHTEA VQAALAKYKE RKMPMPKRR SVLVHSSVET YTPDTSSAS EDEGLRRPG RLTSTPLQSH SSVEPWLDLV IQGSSTSSSA SSTSSHPGGR

image



Background

Alternative Names Disco-interacting protein 2 homolog A (DIP2 homolog A)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human DIP2A Protein
Catalog #: PEH10395

