


Product Name: Recombinant Human UBN1 Protein
Catalog #: PEH10399



Summary

Name	UBN1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NPG3
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LLQAATAMDS LTDLLEHLL SESPEGSPFR DMDDGSDSLG VGLDQEFRQP SSLPEGLPAP LEKRVKELAQ AARAAEGESR QKFFTQDING ILLDIEAQTR ELSSQVRSGV YAYLASFLPC SKDALLKRAR KLHLYEQGGR LKEPLQKLKE AIGRAMPEQM AKYQDECQAH TQAKVAKMLE EEKDKEQRDR ICSDEEEDDEE KGGRRIMGPR KKFQWNDEIR ELLCQVVKIK LESQDLERNN KAQAWEDCVK GFLDAEVKPL

image



Background

Alternative Names	Ubiquitin-1 (HIRA-binding protein) (Protein VT4) (Ubiquitously expressed nuclear protein)
--------------------------	---

Product Name: Recombinant Human UBN1 Protein
Catalog #: PEH10399



Note

For Research Use Only , Not for Diagnostic Use