


Product Name: Recombinant Human NAP1L3 Protein
Catalog #: PEH10405



Summary

Name	NP1L3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99457
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEEENPKENP EVKAEKEVP KEIPEVKDEE KEVPKEIPEV KAEKADSKD CMEATPEVKE DPKEVPQVKA DDKEQPKATE AKARAAVRET HKRVPEERLQ DSVDLKRARK GKPKREDPKG IPDYWLIVLK NVDKLGPMIQ KYDEPILKFL SDVSLKFSKP GQPVSYTFEF HFLPNPYFRN EVLVKTYIIK AKPDHNDPFF SWGWEIEDCK GCKIDWRRGK DVTVTQTQSR TTATGEIEIQ PRVVPNASFF NFFSPPEIPM IGGLEPREDA ILDEDFEIGQ ILHDNVILKS IYYTGEVNG

image



Background

Alternative Names Nucleosome assembly protein 1-like 3

Product Name: Recombinant Human NAP1L3 Protein
Catalog #: PEH10405



Note

For Research Use Only , Not for Diagnostic Use