


Product Name: Recombinant Human CCNT1 Protein
Catalog #: PEH10407



Summary

Name	CCNT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O60563
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GKRWLSSQPS FKLEPTQGHR TSENALTGV DHSLPQDGSN AFISQKQNSK SVPSAKVSLK EYRAKHAEEL AAQKRQLENM EANVKSQYAY AAQNLLSHHD SHSSVILKMP IEGSENPERP FLEKADKTAL KMRIPVAGGD KAASSKPEEI KMRIKVHAAA DKHNSVEDSV TKSREHKEKH KTHPSNHHHH HNHSHKHSH SQLPVGTGNK RPGDPKHSSQ TSNLAHKTY

image



Background

Alternative Names Cyclin-T1 (CycT1) (Cyclin-T)

Note

Product Name: Recombinant Human CCNT1 Protein
Catalog #: PEH10407



For Research Use Only , Not for Diagnostic Use