


Product Name: Recombinant Human INTS8 Protein
Catalog #: PEH10409



Summary

Name	INT8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q75QN2
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VLEAALKLNK DLYVHTMRTL DLLAMEPGMV NGETESSTAG LKVKTEEMQC QVCYDLGAAY FQQGSTNSAV YENAREKFFR TKELIAEIGS LSLHCTIDEK RLAGYCQACD VLPSSDSTS QQLTPYSQVH ICLRSGNYQE VIQIFIEDNL TLSLPVQFRQ SVLRELFKKA QQGNEALDEI CFKVCACNTV RDILEGRTIS VQFNQLFLRP NKEKIDFLE VCSRSVNLEK ASES�KGNMA AFLKNVCLGL

image



Background

Alternative Names Integrator complex subunit 8 (Int8) (Protein kaonashi-1)

Note

Product Name: Recombinant Human INTS8 Protein
Catalog #: PEH10409



For Research Use Only , Not for Diagnostic Use