

Product Name: Recombinant Human SMARCD2 Protein
Catalog #: PEH10410



Summary

Name	SMRD2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92925
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HVISVDPNDQ KKTACYDIDV EVDDPLKAQM SNFLASTTNQ QEIASLDVKI HETIESINQL KTQRDFMLSF STDPQDFIQE WLRSQRRDLK IITDVIGNPE

image



Background

Alternative Names	SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily D member 2 (60 kDa BRG-1/Brm-associated factor subunit B) (BRG1-associated factor 60B) (BAF60B)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use