


Product Name: Recombinant Human PDS5A Protein
Catalog #: PEH10415



Summary

Name	PDS5A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q29RF7
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KNIDPVKNKE INSDQATQGN ISSDRGKKRT VTAAGAENIQ QKTDEKVDDES GPPAPSKPRR GRRPKSESQG NATKNDDLNK PINKGRKRAA VGQESPGGLE AGNAKAPKLQ DLAKKAAPAE

image



Background

Alternative Names	Sister chromatid cohesion protein PDS5 homolog A (Cell proliferation-inducing gene 54 protein) (Sister chromatid cohesion protein 112) (SCC-112)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use