


Product Name: Recombinant Human ZNF408 Protein
Catalog #: PEH10416



Summary

Name	ZN408_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H9D4
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GKKALQLARE PRLGLDLGWN PSGEGCTQGL KDVPPEPTRD ILALKSLPRG LALGPSLAKE QRLGVWCVD PLQPGLLWGP LEEESASKEK GEGVKPRQEE NLSLGPWGDV CACEQSSGWT SLVQRGRLES EGNVAPVRIS ERLHLQVYQL VLPGSELLW PQPSSEGPSL TQPGLDKEAA VAVVTEVES VQQEVASPG

image



Background

Alternative Names Zinc finger protein 408 (PR domain zinc finger protein 17)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZNF408 Protein
Catalog #: PEH10416

