


Product Name: Recombinant Human CWC25 Protein
Catalog #: PEH10417



Summary

Name	CWC25_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NXE8
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SWHPQTLRNV EKVWKAQKH EAERKKIEEL QRELREERAR EEMQRYAEDV GAVKKKEEKL DWMYQGPGGM

image



Background

Alternative Names	Pre-mRNA-splicing factor CWC25 homolog (Coiled-coil domain-containing protein 49) (Spliceosome-associated protein homolog CWC25)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use