


Product Name: Recombinant Human RNF169 Protein
Catalog #: PEH10418




Summary

Name	RN169_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NCN4
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ECARRSQPER CRPRRDGGAA AAGPRPEQEP RAAPAEPDFI FRAPIKLSKP GELREEYESL RKLREEKLQE EKPSDQIHK LLPEDTETGK RKMDEQKKRD EPLVLKTNLE RCPARLSDSE NEEPSRGQMT QTHRSAFVSK NNSYSLAFLA GKLNSKVERS QSCSDTAQER AKSRVRAVPG NKAKVTTMTP ASNPIIGVLL STQNNRCVSA PDLTIEKRLP FSSLSSLASL HKPERSVSPE SNDSISEELN

image



SDS



WB

Background

Alternative Names	E3 ubiquitin-protein ligase RNF169 (EC 2.3.2.27) (RING finger protein 169) (RING-type E3 ubiquitin transferase RNF169)
--------------------------	--

Note

Product Name: Recombinant Human RNF169 Protein
Catalog #: PEH10418



For Research Use Only , Not for Diagnostic Use