


Product Name: Recombinant Human ZNF644 Protein
Catalog #: PEH10419



Summary

Name	ZN644_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H582
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LQLSTPQKAS QHQVLFLLSD VAHAKNPHTS NKKLPTSASV GCDIQNSVGS NIKSDGTLIN QVEVGEDGED LLVKDDCVNT VTGISSGTDG FRSENDTNWD PQKEFIQFLM TNEETVDKAP PHSKIGLEKK RKRKMDVSKI TRYTEDCFSD SNCVPNKSKM QEVDLFLEQNE ELQAVDSQKY

image



Background

Alternative Names Zinc finger protein 644 (Zinc finger motif enhancer-binding protein 2) (Zep-2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZNF644 Protein
Catalog #: PEH10419

