


Product Name: Recombinant Human ZNF382 Protein
Catalog #: PEH10420



Summary

Name	ZN382_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96SR6
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KVLKNISELV IRNISPIKEK FGDSTGWEKS LLNTKHEKIH PAVNLHKQTE RVLSGKQELI QHQKVQAPEQ PFDHNECEKS

image



Background

Alternative Names Zinc finger protein 382 (KRAB/zinc finger suppressor protein 1) (KS1)
(Multiple zinc finger and krueppel-associated box protein KS1)

Note

For Research Use Only , Not for Diagnostic Use