


Product Name: Recombinant Human ZHX3 Protein
Catalog #: PEH10421



Summary

Name	ZHX3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H4I2
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FTPTKYKERA PEQLRALESS FAQNPLPLDE ELDRLRSETK MTRREIDSWF SERRKKVNAE ETKKAEENAS QEEEEAAEDE GGEEDLASEL RVSGENGSL MPSSHILAER KVSPKINLK NLRVTEANGR NEIPGLGACD PEDDESINKLA EQLPGKVSCK KTAQQRHLLR QLFVQTQWPS

image



Background

Alternative Names	Zinc fingers and homeoboxes protein 3 (Triple homeobox protein 1) (Zinc finger and homeodomain protein 3)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZHX3 Protein
Catalog #: PEH10421

