


Product Name: Recombinant Human NFRKB Protein
Catalog #: PEH10424



Summary

Name	NFRKB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6P4R8
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PEEREWRTQQ RYLKVLREVK EECGDTALSS DEEDLSSWLP SSPARSPSPA VPLRVVPTLS TTDMMKTADKV ELGDSDLKIM LKKHHEKRKH

image



Background

Alternative Names Nuclear factor related to kappa-B-binding protein (DNA-binding protein R kappa-B) (INO80 complex subunit G)

Note

For Research Use Only , Not for Diagnostic Use