


Product Name: Recombinant Human SKOR1 Protein
Catalog #: PEH10430



Summary

Name	SKOR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P84550
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PCGPLGGPAP AKVFAPERDE HVKSAVALG PAASYVCTPE AHEPDKEDNH SPADDLETRK SYPDQRSISQ PSPANTDRGE DGLTLDVTGT HLVEKDIENTL AREELQKLL EQMELRKKLE REFQSLKDNF QDQMKRELAY REEMVQQLQI VRDTLCNELD QERKARYAIQ QKLKEAHDAL

image



Background

Alternative Names	SKI family transcriptional corepressor 1 (Functional Smad-suppressing element on chromosome 15) (Fussel-15) (LBX1 corepressor 1) (Ladybird homeobox corepressor 1)
--------------------------	--

Note

Product Name: Recombinant Human SKOR1 Protein
Catalog #: PEH10430



For Research Use Only , Not for Diagnostic Use