


Product Name: Recombinant Human PYGO1 Protein
Catalog #: PEH10433



Summary

Name	PYGO1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y3Y4
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YNNALSQNVN MPNQHFRQNP AENFSQIPPQ NASQVSNPDL ASNFVPGNNS NFTSPLESNH SFIPPNTFG QAKAPPPKQD FTQGATKNTN QNSSAHPPHL NMDDTVNQSN IELKNVNRNN AVNQENSRSS STEATNNNPA NGTQNKPRQP RGAADACTTE KSNKSSLHPN RHGHSSSDPV YPCGICTNEV NDDQDAILCE ASCQKWFHRI

image



Background

Alternative Names Pygopus homolog 1

Note

Product Name: Recombinant Human PYGO1 Protein
Catalog #: PEH10433



For Research Use Only , Not for Diagnostic Use