


Product Name: Recombinant Human RCBTB1 Protein
Catalog #: PEH10434



Summary

Name	RCBT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NDN9
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DEVFVFLNY SNCLGTGDNQ STLVPPKLEG LCGKKIKSLS YGSGPHVLLS TEDGVVYAWG HNGYSQLGNG TTNQGIAPVQ VCTNLLIKQV VEVACGSHHS MALAADGEVF AWGYNNCGQV GSGSTANQPT PRKVTNCLHI

image



Background

Alternative Names	RCC1 and BTB domain-containing protein 1 (Chronic lymphocytic leukemia deletion region gene 7 protein) (CLL deletion region gene 7 protein) (Regulator of chromosome condensation and BTB domain-containing protein 1)
--------------------------	--

Note

Product Name: Recombinant Human RCBTB1 Protein
Catalog #: PEH10434



For Research Use Only , Not for Diagnostic Use