

Product Name: Recombinant Human SNRNP27 Protein
Catalog #: PEH10441



Summary

| | |
|------------------------------|--|
| Name | SNR27_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q8WVK2 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 11 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | RDEEKKETKE TSKERQITE EDLEGKTEEE IEMMKLMGFA SFDSTKGKKV DGSVNAYAIN VSQKRKYRQY MNRKGGFNRP |

image



Background

| | |
|--------------------------|---|
| Alternative Names | U4/U6.U5 small nuclear ribonucleoprotein 27 kDa protein (U4/U6.U5 snRNP 27 kDa protein) (U4/U6.U5-27K) (Nucleic acid-binding protein RY-1) (U4/U6.U5 tri-snRNP-associated 27 kDa protein) (27K) (U4/U6.U5 tri-snRNP-associated protein 3) |
|--------------------------|---|

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SNRNP27 Protein
Catalog #: PEH10441

