


Product Name: Recombinant Human DBX2 Protein
Catalog #: PEH10443



Summary

Name	DBX2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZNG2
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LEKMFQKQKY ISKTDRKKLA INLGLKESQV KIWFQNRRMK WRNSKEKEVL SNRCIQEVGL QEDPLRSAL GFPSPCPSIW DVPQQHSSPR WRENSPEPSE RLIQESSGAP PPEANSLQGA LYLCSEEEAG

image



Background

Alternative Names Homeobox protein DBX2 (Developing brain homeobox protein 2)

Note

For Research Use Only , Not for Diagnostic Use