


Product Name: Recombinant Human DDX47 Protein
Catalog #: PEH10447



Summary

Name	DDX47_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H0S4
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YYIFIPSKFK DTYLVYILNE LAGNSFMIFC STCNNTQRTA LLLRNLGFTA IPLHGQMSQS KRLGSLNKFK AKARSILLAT DVASRGLDIP HVDVVVNFDI PTHSKDYIHR VGR TARAGRS GKAITFVTQY DVELFQRIEH LIGKKLPGFP TQDDEVMMILT ERVAEAQRFA RMELREHGEK KKRSREDAGD

image



Background

Alternative Names	Probable ATP-dependent RNA helicase DDX47 (EC 3.6.4.13) (DEAD box protein 47)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human DDX47 Protein
Catalog #: PEH10447

