



Product Name: Recombinant Human BRCA1 Protein
Catalog #: PEH10448



Summary

Name	BRCA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P38398
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KEWLNKQKLPC SENPRDTEVD PWITLNSSIQ KVNEWFSRSD ELLGSDDSHD GESENAKVA DVLVDLNEVD EYSGSSEKID LLASDPHEAL ICKSERVHSK SVESNIEDKI FGKTYRKKAS LPNLSHVTEN LIIGAFVTEP QIIQERPLTN KLKRKRRTS GLHPEDFIKK ADLAVQKTPE MINQGTNQTE QNGQVMNITN

image



Background

Alternative Names Breast cancer type 1 susceptibility protein (EC 2.3.2.27) (RING finger protein 53) (RING-type E3 ubiquitin transferase BRCA1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human BRCA1 Protein
Catalog #: PEH10448

