


Product Name: Recombinant Human ATM Protein
Catalog #: PEH10449



Summary

Name	ATM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13315
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DTNGNLMEVE DQSSMNLFND YPDSSVSDAN EPGESQSTIG AINPLAEEYL SKQDLLFLDM LKFLCLCVTT AQTNTVSFRA ADIRRKLLML IDSSTLEPTK

image



Background

Alternative Names Serine-protein kinase ATM (EC 2.7.11.1) (Ataxia telangiectasia mutated) (A-T mutated)

Note

For Research Use Only , Not for Diagnostic Use