


Product Name: Recombinant Human PPA1 Protein
Catalog #: PEH10455



Summary

Name	IPYR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15181
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PFHDIPIYAD KDV FHMVVEV PRWSNAKMEI ATKDPLNPIK QDVKKGKLRY VANLFPYKGY IWN YGAIPQT WEDPGHNDKH TGCCGDNDPI DVCEIGSKVC ARGEIIGVKV LGILAMIDEG ETDWKVIIN VDDPDAANYN DINDVKRLKP GYLEATVDWF RRYKVPDGKP ENEFAFNAEF KDKDFAIDII KSTHDPHWKAL VTKKTNGKGI SCMNTTLES PFKCDPDAAR AIVDALPPPC

image



Background

Alternative Names	Inorganic pyrophosphatase (EC 3.6.1.1) (Pyrophosphate phospho-hydrolase) (PPase)
--------------------------	--

Note

Product Name: Recombinant Human PPA1 Protein
Catalog #: PEH10455



For Research Use Only , Not for Diagnostic Use